

Phospho-HSP90AB1 (S226) Antibody



PACO02822

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

HSP90AB1

Uniprot

P08238

Synonyms:

HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84

Immunogen:

Synthesized peptide derived from human HSP90Beta around the phosphorylation site of S226.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western blot analysis of VEC using p-HSP90Beta (S226) antibody.

