

OR2A1; Antibody



PACO02837

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2A1

Uniprot

Q8NGT9

Synonyms:

OR2A1; OR2A42; Olfactory receptor 2A1/2A42; Olfactory receptor OR7-16; Olfactory receptor OR7-19

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2A42.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COLO cells using Olfactory receptor 2A42 Polyclonal Antibody.

