

Phospho-IL2RB (Y364) Antibody



PACO02845

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

IL2RB

Uniprot

P14784

Synonyms:

IL2RB; Interleukin-2 receptor subunit beta; IL-2 receptor subunit beta; IL-2R subunit beta; IL-2RB; High affinity IL-2 receptor subunit beta; p70-75; p75; CD antigen CD122

Immunogen:

Synthesized peptide derived from human IL-2RBeta around the phosphorylation site of Y364.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293T cells using Phospho-IL-2RBeta (Y364) Polyclonal Antibody.

