

OR52K1 Antibody



PACO02859

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR52K1

Uniprot

Q8NGK4

Synonyms:

OR52K1; Olfactory receptor 52K1; Olfactory receptor OR11-8

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 52K1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293T cells using Olfactory receptor 52K1 Polyclonal Antibody.

