

OR10V1 Antibody



PACO02863

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR10V1

Uniprot

Q8NGI7

Synonyms:

OR10V1; Olfactory receptor 10V1; Olfactory receptor OR11-256


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 10V1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Olfactory receptor 10V1 Polyclonal Antibody.