

LAMA5 Antibody



PACO02865

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

LAMA5

Uniprot

O15230

Synonyms:

LAMA5; KIAA0533; KIAA1907; Laminin subunit alpha-5; Laminin-10 subunit alpha; Laminin-11 subunit alpha; Laminin-15 subunit alpha


Immunogen:

Synthesized peptide derived from the Internal region of human Laminin alpha-5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAW HepG2 A549 cells using Laminin alpha-5 Polyclonal Antibody.