

Phospho-LCP2 (Y128) Antibody



PACO02868

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

LCP2

Uniprot

Q13094

Synonyms:

LCP2; Lymphocyte cytosolic protein 2; SH2 domain-containing leukocyte protein of 76 kDa; SLP-76 tyrosine phosphoprotein; SLP76

Immunogen:

Synthesized peptide derived from human SLP-76 around the phosphorylation site of Y128.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of K562 cells using Phospho-SLP-76 (Y128) Polyclonal Antibody.

