

LIMK1 Antibody



PACO02874

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

LIMK1

Uniprot

P53667

Synonyms:

LIMK1; LIMK; LIM domain kinase 1; LIMK-1


Immunogen:

Synthesized peptide derived from human LIMK-1 around the non-phosphorylation site of T508.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using LIMK-1 Polyclonal Antibody.