

LIMK2 Antibody



PACO02876

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

LIMK2

Uniprot

P53671

Synonyms:

LIMK2; LIM domain kinase 2; LIMK-2


Immunogen:

Synthesized peptide derived from human LIMK-2 around the non-phosphorylation site of S283.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 3T3 cells using LIMK-2 Polyclonal Antibody.