

PACO02877

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

LIPC

Uniprot

P11150

Synonyms:

LIPC; HTGL; Hepatic triacylglycerol lipase; HL; Hepatic lipase; Lipase member C


Immunogen:

Synthesized peptide derived from the Internal region of human Hepatic Lipase.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Hepatic Lipase Polyclonal Antibody.