

Phospho-MCL1 (S159) Antibody



PACO02894

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

MCL1

Uniprot

Q07820

Synonyms:

MCL1; BCL2L3; Induced myeloid leukemia cell differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl2-L-3; Bcl-2-related protein EAT/mcl1; mcl1/EAT

Immunogen:

Synthesized peptide derived from human Mcl-1 around the phosphorylation site of S159.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western blot analysis of K562 using p-Mcl-1 (S159) antibody.

