

MNDA Antibody



PACO02908

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

MNDA

Uniprot

P41218

Synonyms:

MNDA; Myeloid cell nuclear differentiation antigen

Immunogen:

Synthesized peptide derived from the C-terminal region of human MNDA.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 3T3 cells using MNDA Polyclonal Antibody.

