

PACO02912

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MSN/RDX/EZR

Uniprot

P26038/P35241/P15311

Synonyms:

MSN; Moesin; Membrane-organizing extension spike protein; RDX; Radixin; EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81


Immunogen:

Synthesized peptide derived from human Moesin/Ezrin/Radixin around the non-phosphorylation site of T558.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 3T3 cells using Moesin/Ezrin/Radixin Polyclonal Antibody.