

MYB Antibody



PACO02917

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

MYB

Uniprot

P10242

Synonyms:

MYB; Transcriptional activator Myb; Proto-oncogene c-Myb


Immunogen:

Synthesized peptide derived from human c-Myb around the non-phosphorylation site of S12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of KB K562 cells using c-Myb Polyclonal Antibody.