

MYBL2 Antibody



PACO02919

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

MYBL2

Uniprot

P10244

Synonyms:

MYBL2; BMYB; Myb-related protein B; B-Myb; Myb-like protein 2

Immunogen:

Synthesized peptide derived from human B-Myb around the non-phosphorylation site of S577.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 3T3 cells using B-Myb Polyclonal Antibody.

