

PABPC3 Antibody



PACO02946

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IHC:1:200-1:1000

Protein Background:

Gene ID:

PABPC3

Uniprot

Q9H361

Synonyms:

PABPC3; PABP3; PABPL3; Polyadenylate-binding protein 3; PABP-3; Poly(A)-binding protein 3; Testis-specific poly(A)-binding protein

Immunogen:

Synthesized peptide derived from the Internal region of human PABP3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 3T3 cells using PABP3 Polyclonal Antibody.

