

PAK1/PAK2/PAK3 Antibody



PACO02950

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PAK1/PAK2/PAK3

Uniprot

Q13153/Q13177/O75914

Synonyms:

PAK1; Serine/threonine-protein kinase PAK 1; Alpha-PAK; p21-activated kinase 1; PAK-1; p65-PAK; PAK2; Serine/threonine-protein kinase PAK 2; Gamma-PAK; PAK65; S6/H4 kinase; p21-activated kinase 2; PAK-2; p58; PAK3; OPHN3; Serine/threonine-p

Immunogen:

Synthesized peptide derived from human PAKalpha/Beta/gamma around the non-phosphorylation site of T423/402/421.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 3T3 cells using PAKalpha/Beta/gamma Polyclonal Antibody.

