

Phospho-PARK2 (S131) Antibody



PACO02954

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

PRKN

Uniprot

O60260

Synonyms:

PARK2; PRKN; E3 ubiquitin-protein ligase parkin; Parkinson juvenile disease protein 2; Parkinson disease protein 2


Immunogen:

Synthesized peptide derived from human Parkin around the phosphorylation site of S131.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of 293T AD293 HeLa using p-Parkin (S131) antibody.