

# PIPOX Antibody



PACO02961

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

PIPOX

**Uniprot**

Q9P0Z9

**Synonyms:**

PIPOX; LPIPOX; PSO; Peroxisomal sarcosine oxidase; PSO; L-pipecolate oxidase; L-pipecolic acid oxidase

**Immunogen:**

Synthesized peptide derived from the Internal region of human PIPOX.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 3T3 cells using PIPOX Polyclonal Antibody.

