

PLA1A Antibody



PACO02962

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PLA1A

Uniprot

Q53H76

Synonyms:

PLA1A; NMD; PSPLA1; Phospholipase A1 member A; Phosphatidylserine-specific phospholipase A1; PS-PLA1

Immunogen:

Synthesized peptide derived from the C-terminal region of human PLA1A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 3T3 cells using PLA1A Polyclonal Antibody.

