

PEG3 Antibody



PACO02966

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

PEG3

Uniprot

Q9GZU2

Synonyms:

PEG3; KIAA0287; ZSCAN24; Paternally-expressed gene 3 protein; Zinc finger and SCAN domain-containing protein 24

Immunogen:

Synthesized peptide derived from the Internal region of human PEG3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of A549 cells using PEG3 Polyclonal Antibody.

