

# PGA3/PGA4/PGA5 Antibody



PACO02968

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

PGA3/PGA4/PGA5

**Uniprot**

P00790

**Synonyms:**

Pepsin A; PGA3; PGA4; PGA5

**Immunogen:**

Synthesized peptide derived from the Internal region of human Pepsin A.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of K562 cells using Pepsin A Polyclonal Antibody.

