

PGD Antibody



PACO02969

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

PGD

Uniprot

P52209

Synonyms:

PGD; PGDH; 6-phosphogluconate dehydrogenase; decarboxylating

Immunogen:

Synthesized peptide derived from the Internal region of human PGD.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of A549 cells using PGD Polyclonal Antibody.

