

Phospho-PGR (S294) Antibody



PACO02970

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

PGR

Uniprot

P06401

Synonyms:

PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3


Immunogen:

Synthesized peptide derived from human PR around the phosphorylation site of S294.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Phospho-PR (S294) Polyclonal Antibody.