

Phospho-PLA2G4A (S505) Antibody



PACO02974

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

PLA2G4A

Uniprot

P47712

Synonyms:

PLA2G4A; CPLA2; PLA2G4; Cytosolic phospholipase A2; cPLA2; Phospholipase A2 group IVA

Immunogen:

Synthesized peptide derived from human cPLA2 around the phosphorylation site of S505.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western blot analysis of VEC using p-cPLA2 (S505) antibody.

