

Phospho-PLK1 (T210) Antibody



PACO02976

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

PLK1

Uniprot

P53350

Synonyms:

PLK1; PLK; Serine/threonine-protein kinase PLK1; Polo-like kinase 1; PLK-1; Serine/threonine-protein kinase 13; STPK13


Immunogen:

Synthesized peptide derived from human Plk around the phosphorylation site of T210.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of 293 using p-Plk (T210) antibody.