

Phospho-PPP2CA (Y307) Antibody



PACO02992

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

PPP2CA

Uniprot

P67775

Synonyms:

PPP2CA; Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform; PP2A-alpha; Replication protein C; RP-C

Immunogen:

Synthesized peptide derived from human PP2A-Calpha around the phosphorylation site of Y307.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western blot analysis of SKOV3 using p-PP2A-Calpha (Y307) antibody.

