

Phospho-PRKD1 (S910) Antibody



PACO03009

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

PRKD1

Uniprot

Q15139

Synonyms:

PRKD1; PKD; PKD1; PRKCM; Serine/threonine-protein kinase D1; Protein kinase C mu type; Protein kinase D; nPKC-D1; nPKC-mu

Immunogen:

Synthesized peptide derived from human PKD1 around the phosphorylation site of S910.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of A431 cells using Phospho-PKD1 (S910) Polyclonal Antibody.

