

PAK7/PAK6 Antibody



PACO03032

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

PAK5/PAK6

Uniprot

Q9P286/Q9NQ5

Synonyms:

PAK7; KIAA1264; PAK5; Serine/threonine-protein kinase PAK 7; p21-activated kinase 5; PAK-5; p21-activated kinase 7; PAK-7; PAK6; PAK5; Serine/threonine-protein kinase PAK 6; PAK-5; p21-activated kinase 6; PAK-6


Immunogen:

Synthesized peptide derived from human PAK5/6 around the non-phosphorylation site of S602/S560.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO cells using PAK5/6 Polyclonal Antibody.