

# Phospho-REL (S503) Antibody



PACO03063

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

REL

**Uniprot**

Q04864

**Synonyms:**

REL; Proto-oncogene c-Rel

**Immunogen:**

Synthesized peptide derived from human c-Rel around the phosphorylation site of S503.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A