

# SDC4 Antibody



PACO03075

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

SDC4

**Uniprot**

P31431

**Synonyms:**

SDC4; Syndecan-4; SYND4; Amphiglycan; Ryudocan core protein

**Immunogen:**

Synthesized peptide derived from human Syndecan-4 around the non-phosphorylation site of S179.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A