

# VAMP1/VAMP2/VAMP3 Antibody



PACO03098

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

VAMP1/VAMP2/VAMP3

**Uniprot**

P23763/P63027/Q15836

**Synonyms:**

VAMP1; SYB1; Vesicle-associated membrane protein 1; VAMP-1; Synaptobrevin-1; VAMP2; SYB2; Vesicle-associated membrane protein 2; VAMP-2; Synaptobrevin-2; VAMP3; SYB3; Vesicle-associated membrane protein 3; VAMP-3; Cellubrevin; CEB; Synaptob

**Immunogen:**

Synthesized peptide derived from the Internal region of human VAMP-1/2/3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A