

Phospho-TGFBR1 (S165) Antibody



PACO03108

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:50-1:300

Protein Background:**Gene ID:**

TGFBR1

Uniprot

P36897

Synonyms:

TGFBR1; ALK5; SKR4; TGF-beta receptor type-1; TGFR-1; Activin A receptor type II-like protein kinase of 53kD; Activin receptor-like kinase 5; ALK-5; ALK5; Serine/threonine-protein kinase receptor R4; SKR4; TGF-beta type I receptor; Transfor

Immunogen:

The antiserum was produced against Synthesized peptide derived from human TGF beta Receptor I around the phosphorylation site of Ser165. AA range:131-180.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A