

Phospho-TH (S62) Antibody



PACO03110

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TH

Uniprot

P07101

Synonyms:

TH; TYH; Tyrosine 3-monooxygenase; Tyrosine 3-hydroxylase; TH

Immunogen:

Synthesized peptide derived from human TH around the phosphorylation site of S62.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A