

# Phospho-VIM (S56) Antibody



PACO03127

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

VIM

**Uniprot**

P08670

**Synonyms:**

VIM; Vimentin

**Immunogen:**

Synthesized peptide derived from human Vimentin around the phosphorylation site of S56.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A