

# Phospho-WEE1 (S642) Antibody



PACO03128

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

WEE1

**Uniprot**

P30291

**Synonyms:**

WEE1; Wee1-like protein kinase; WEE1hu; Wee1A kinase

**Immunogen:**

Synthesized peptide derived from human Wee1 around the phosphorylation site of S642.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A