

Phospho-Histone H2A (T121) Antibody



PACO03143

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

HIST1H2AG/HIST1H2AB/HIST1H2AD/HIST1H2AH/HIST2H2AA3/HIST3H2A

Uniprot

P0C0S8/P04908/P20671/Q96KK5/Q6FI13/Q7L7L0

Synonyms:

HIST1H2AG; H2AFP; HIST1H2AI; H2AFC; HIST1H2AK; H2AFD; HIST1H2AL; H2AFI; HIST1H2AM; H2AFN; Histone H2A type 1; H2A.1; Histone H2A/p; HIST1H2AB; H2AFM; HIST1H2AE; H2AFA; Histone H2A type 1-B/E; Histone H2A.2; Histone H2A/a; Histone H2A/m; HIS

Immunogen:

Synthesized peptide derived from human Histone H2A around the phosphorylation site of T121.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A