

Phospho-CASP6 (S257) Antibody



PACO03150

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

CASP6

Uniprot

P55212

Synonyms:

CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2

Immunogen:

Synthesized peptide derived from human Caspase-6 around the phosphorylation site of S257.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A