

TECR Antibody



PACO03182

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

TECR

Uniprot

Q9NZ01

Synonyms:

TECR; GPSN2; SC2; Trans-2; 3-enoyl-CoA reductase; TER; Synaptic glycoprotein SC2

Immunogen:

Synthesized peptide derived from the C-terminal region of human GPSN2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A