

Phospho-PGC (Ser15) Antibody



PACO03212

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Plants	PGC
Source:	Uniprot
Rabbit	P20142
Isotype:	Synonyms:
IgG	PGC; Gastricsin; Pepsinogen C
Applications:	Immunogen:
WB	Synthetic peptide.
Recommended dilutions:	Storage:
	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A