

TLR13 Antibody



PACO03217

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:200-1:500

Protein Background:

Gene ID:

Tlr13

Uniprot

Q6R5N8

Synonyms:

Toll-like receptor 13

Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A