

RPS6KB1 Antibody



PACO03225

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RPS6KB1

Uniprot

P23443

Synonyms:

Ribosomal protein S6 kinase beta-1

Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A