

HAO1 Antibody



PACO03226

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HAO1

Uniprot

Q9UJM8

Synonyms:

Hydroxyacid oxidase 1

Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A