

H2AFZ Antibody



PACO03231

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

Protein Background:

Gene ID:

H2AFZ

Uniprot

P0C0S5

Synonyms:

H2AFZ; H2AZ; Histone H2A. Z; H2A/z

Immunogen:

Synthetic peptide.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A