

# IRAK1 Antibody



PACO03241

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

IRAK1

**Uniprot**

P51617

**Synonyms:**

IRAK1; IRAK; Interleukin-1 receptor-associated kinase 1; IRAK-1

**Immunogen:**

Synthesized peptide derived from human IRAK-1 around the non-phosphorylation site of T209.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A