

OR51B2 Antibody



PACO03277

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR51B2

Uniprot

Q9Y5P1

Synonyms:

OR51B2; OR51B1P; Olfactory receptor 51B2; Odorant receptor HOR5'beta3; Olfactory receptor 51B1

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 51B2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A