

HSPA7 Antibody



PACO03287

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

HSPA7

Uniprot

P48741

Synonyms:

HSPA7; HSP70B; Putative heat shock 70 kDa protein 7; Heat shock 70 kDa protein B; HSPA6; HSP70B'; Heat shock 70 kDa protein 6; Heat shock 70 kDa protein B'

Immunogen:

Synthesized peptide derived from the Internal region of human HSP77/76.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A