

# ADAR Antibody



PACO03335

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

ADAR

**Uniprot**

P55265

**Synonyms:**

ADAR; ADAR1; DSRAD; G1P1; IFI4; Double-stranded RNA-specific adenosine deaminase; DRADA; 136 kDa double-stranded RNA-binding protein; p136; Interferon-inducible protein 4; IFI-4; K88DSRBP

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human ADAR1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A