

CIB2 Antibody



PACO03346

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

CIB2

Uniprot

O75838

Synonyms:

CIB2; KIP2; Calcium and integrin-binding family member 2; Kinase-interacting protein 2; KIP 2

Immunogen:

Synthesized peptide derived from the Internal region of human CIB2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A