

# Phospho-CENPA (S7) Antibody



PACO03354

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

CENPA

**Uniprot**

P49450

**Synonyms:**

CENPA; Histone H3-like centromeric protein A; Centromere autoantigen A; Centromere protein A; CENP-A

**Immunogen:**

Synthesized peptide derived from human CENP-A around the phosphorylation site of S7.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A